



MANUAL
WHEELCHAIRS
OPTIONS AND
ACCESSORIES
GUIDE



TABLE OF CONTENTS



Seat upholstery	4	Spokeguards.....	53
Backrest	7	Wheel locks	59
Push handles	11	Accessories	62
Back upholstery	13	Colour options	66
Sideguards and armrests	14	Useful resources.....	72
Legrests and footrests	19		
Frame versions	28		
Centre of gravity.....	32		
Frame accessories	36		
Forks.....	41		
Castor wheels.....	43		
Rear wheels	45		
Handrims	50		

SEAT UPHOLSTERY

Seat sling standard - with one utility bag

Standard nylon sling with a utility bag underneath - useful for storing belongings and tools.

Available on:

- | | |
|-------------------------|-----------|
| ✓ Krypton F | KRF020002 |
| ✓ Helium | HUM020002 |
| ✓ Xenon ² FF | XFF020002 |
| ✓ Xenon ² H | XFF020002 |
| ✓ Xenon ² SA | XES020002 |
| ✓ Argon ² | AGN020002 |
| ✓ Neon ² | NSW020002 |
| ✓ Life F | LFE020014 |
| ✓ Life R | LRF020002 |



Similar options

Seat sling standard - with two utility bags

The same style as the above with an added bag.

Available on:

- | | |
|-------------------------|-----------|
| ✓ Krypton F | KRF020003 |
| ✓ Helium | HUM020003 |
| ✓ Xenon ² FF | XFF020003 |
| ✓ Xenon ² H | XFF020003 |
| ✓ Xenon ² SA | XES020003 |
| ✓ Argon ² | AGN020003 |
| ✓ Neon ² | NSW020003 |

Seat sling standard - with one utility bag and catheter bag

The same style as the above with an added catheter bag.

Available on:

- | | |
|-------------------------|-----------|
| ✓ Krypton F | KRF020004 |
| ✓ Helium | HUM020004 |
| ✓ Xenon ² FF | XFF020004 |
| ✓ Xenon ² H | XFF020004 |
| ✓ Xenon ² SA | XES020004 |
| ✓ Argon ² | AGN020004 |
| ✓ Neon ² | NSW020004 |



Seat sling - Krypton R

The new krypton R offers a tension adjustable seat sling that wraps around the frame. Perfect to help achieve the optimum seated position and correct foot position.



Seat sling standard - with extra padding on seat

Additional padding is useful for users that may only be using the chair for very short periods and do not require a cushion.

Available on:

- | | |
|-------------|-----------|
| ✓ Krypton R | KRR020012 |
|-------------|-----------|

Available on:

- | | |
|---------------------|-----------|
| ✓ Neon ² | NSW020050 |
|---------------------|-----------|



Seat sling clip style - airflow (user weight limit 100kg)

The clip style seat sling is our lightest option for the seat upholstery. The straps can be re-tensioned or contoured to suit. The individual straps also promote airflow and ventilation.

Available on:

✓ Helium HUM020010



Similar options

Seat sling clip style - airflow (user weight limit 100kg) - with one utility bag

The same style as the above with a utility bag under the front section of the chair.

Available on:

✓ Helium HUM020005

Seat sling clip style - airflow (user weight limit 100kg) - with two utility bags

The same style as the above with two utility bags under the front section of the chair.

Available on:

✓ Helium HUM020006

Seat sling clip style - airflow (user weight limit 100kg) - with one utility bag and catheter bag

The same style as the above with a utility bag and catheter bag under the front section of the chair.

Available on:

✓ Helium HUM020007

SEAT UPHOLSTERY

Positioning belt with lock

This belt can be used to hold the client in position and prevent them moving forwards in the chair.

Available on:

✓ Krypton R	KRR090018
✓ Krypton F	KRF090018
✓ Helium	HUM090018
✓ Xenon ² FF	XFF090018
✓ Xenon ² H	XFF090018
✓ Xenon ² SA	XES090018
✓ Argon ²	AGN090018
✓ Neon ²	NSW090018
✓ Life F	LFE090018
✓ Life R	LRF090018



Backrest

Backrest angle

Where a backrest angle is an adjustment or option to select on the order form it is listed as a positive or negative angle and is in relation to the the seat frame and not the floor.

This is illustrated in this diagram.

Negative angle squeezes the end user into the chair, which can be useful when the chair is set up with bucket.

Positive angle opens the backrest up and can be used to accommodate those with open hips.



BACKREST

Fixed height backrest (non height adjustable, carbon fibre)

This backrest is not height adjustable but is our lightest offering. This should only be chosen when the user is established and unlikely to change their cushion.

Available on:

- | | |
|-------------------------|-----------|
| ✓ Krypton R | KRR030006 |
| ✓ Krypton F | KRF030007 |
| ✓ Xenon ² FF | XFF030002 |
| ✓ Xenon ² H | XFF030002 |
| ✓ Xenon ² SA | XES030002 |



Angle adjustable back

Back angle is adjustable forwards and backwards, (range is dependant on the chair (see order form)) giving the ability to squeeze the client when a deep bucket is used or open the backrest for those with open hips. (Image shows Xenon)

Available on:

- | | |
|-------------------------|-----------|
| ✓ Krypton R | KRR030012 |
| ✓ Krypton F | FRK030012 |
| ✓ Helium | HUM030012 |
| ✓ Xenon ² FF | XFF030012 |
| ✓ Xenon ² H | XFF030012 |
| ✓ Xenon ² SA | XES030012 |
| ✓ Argon ² | AGN030012 |
| ✓ Neon ² | NSW030012 |
| ✓ Life F | LFE030012 |
| ✓ Life R | LRF030015 |



BACKREST

Back angle adjustable foldable*

The back angle is adjustable forwards and backwards although the range is dependent on the chair (listed on the order form). This gives the ability to squeeze the client when a deep bucket is used or open the backrest for those with open hips. The backrest is also able to fold for transfers and storage.

Available on:

- ✓ Krypton R KRR030015
- ✓ Helium HUM030015
- ✓ Argon² AGN030015

* Rigid chairs only



Double locking (angle adjustable back)

The backrest will lock in the folded position, this means that the end user will be able to grab the support bar at the back to assist transfers. (Image shows Helium)

Available on:

- ✓ Krypton R KRR030102
- ✓ Helium HUM030102
- ✓ Argon² AGN030102
- ✓ Life R LRF030102



Back width one size smaller (SW -2cm)

The backrest tubes sit inside the seat tubes allowing for a tapered backrest. This is popular for women who have slightly wider hips but would like a tighter more supportive fit in the back.

Available on:

- ✓ Helium HUM030103



Standard: non angle adjustable (90°), non folding backrest

This backrest is fixed and will not fold, often chosen when these options are not required and more rigidity is preferred in the backrest. This will be stronger for someone that extends when propelling or who is likely to spasm. It is also lighter as there are less components.

Available on:

- ✓ Xenon² FF XFF030010
- ✓ Xenon² H XFF030010
- ✓ Xenon² SA XES030010
- ✓ Neon² NSW030010
- ✓ Life F LFE030010
- ✓ Life R LRF030001

Half folding back (to front)

The half folding back allows the user or carers to fold down the backrest on a folding chair to reduce the overall package size for stowage.

Available on:

- ✓ Krypton F KRF030021
- ✓ Xenon² FF XFF030031
- ✓ Xenon² H XFF030031
- ✓ Xenon² SA XES030031
- ✓ Neon² NSW030031
- ✓ Life F LFE030031



Folding stabiliser bar

The stabiliser bar may be a good option on wider chairs to offer more rigidity for larger users. Folds up automatically when the chair is folded.

Available on:

- ✓ Krypton F KRF030402
- ✓ Xenon² FF XFF030402
- ✓ Xenon² H XFF030402
- ✓ Xenon² SA XES030402
- ✓ Neon² NSW030402
- ✓ Life F LFE030402



BACKREST

Stabilising Bar Swing Away

The stabiliser bar may be a good option on wider chairs to offer more rigidity for larger users. Swings away when the chair is folded.

Available on:

- ✓ Xenon² SA XES030401
- ✓ Neon² NSW030401

Oval shaped backrest (non height adjustable)

This ergonomically shaped bar is the perfect grab point for transfers into a car. It also has a sleeker look than the standard support bar.

Available on:

- ✓ Krypton R KRR030006
- ✓ Helium HUM030003



Similar options

Stabiliser bar - close setting (9.2cm)

The oval shaped backrest bar can be set at either 9.2cm or 12.4cm (opposite) away from the backrest. If the end user is likely to be submerged into the tension adjustable upholstery or if they are using a hard back then the 12.4cm wide option will offer more space to accommodate for this.

Available on:

- ✓ Helium HUM030403

Stabiliser bar - wide setting (12.4cm)

Available on:

- ✓ Helium HUM030404

PUSH HANDLES

Long push handles

These push handles are great for someone that is being pushed the majority of the time.

Available on:

- | | |
|-------------------------|-----------|
| ✓ Krypton R | KRR030201 |
| ✓ Krypton F | KRF030201 |
| ✓ Helium | HUM030201 |
| ✓ Xenon ² FF | XFF030201 |
| ✓ Xenon ² H | XFF030201 |
| ✓ Xenon ² SA | XES030201 |
| ✓ Argon ² | AGN030201 |
| ✓ Neon ² | NSW030201 |
| ✓ Life F | LFE030201 |
| ✓ Life R | LRF030201 |



Fold-down push handles

The fold-down push handles may be used by pressing the button inside and are perfect for users that do not always require assistance. They are very discreet when folded.

Available on:

- | | |
|-------------------------|-----------|
| ✓ Krypton R | KRR030203 |
| ✓ Krypton F | KRF030203 |
| ✓ Helium | HUM030203 |
| ✓ Xenon ² FF | XFF030203 |
| ✓ Xenon ² H | XFF030203 |
| ✓ Xenon ² SA | XES030203 |
| ✓ Argon ² | AGN030203 |
| ✓ Neon ² | NSW030203 |
| ✓ Life F | LFE030203 |
| ✓ Life R | LRF030203 |




PUSH HANDLES

Height adjustable push handles

Height adjustable push handles are good for users with a low back or chair, or for those who might have multiple carers. They can be dropped down and tucked away when not in use.

Available on:

✓ Krypton R	KRR030204
✓ Krypton F	KRF030204
✓ Helium	HUM030204
✓ Xenon ² FF	XFF030204
✓ Xenon ² H	XFF030204
✓ Xenon ² SA	XES030204
✓ Argon ²	AGN030204
✓ Neon ²	NSW030204
✓ Life F	LFE030204
✓ Life R	LRF030204




Similar options

Omit push handles

No push handles will be provided on the chair.

Available on:

✓ Krypton R	KRR030200
✓ Krypton F	KRF030200
✓ Helium	HUM030200
✓ Xenon	XFF030200
✓ Argon ²	AGN030200
✓ Neon ²	NSW030200



J3 hardware - height adjustable push handles

If a JAY 3 back is going to be mounted and push handles are required then these allow more options around placement and height.

Available on:

✓ Krypton R	KRR030205
✓ Krypton F	KRF030205
✓ Helium	HUM030205
✓ Xenon ² FF	XFF030205
✓ Xenon ² H	XFF030205
✓ Xenon ² SA	XES030205
✓ Argon ²	AGN030205
✓ Neon ²	NSW030205
✓ Life F	LFE030205
✓ Life R	LRF030205

BACK UPHOLSTERY

Exo Upholstery

The Exo is designed for both comfort and breathability. The ribbed curves allow for air to flow and the compressed EVA zote fabric ensures that it is both lightweight and comfortable.

Available on:

- ✓ Krypton R KRR030316
- ✓ Krypton F KRF030316
- ✓ Helium HUM030316
- ✓ Xenon² FF XFF030316
- ✓ Xenon² H XFF030316
- ✓ Xenon² SA XES030316
- ✓ Argon² AGN030316
- ✓ Neon² NSW030316
- ✓ Life F LFE030316
- ✓ Life R LRF030316




Exo Pro Upholstery

The Exo Pro combines the features of the Exo upholstery and more. With 3DX spacer fabric technology from the JAY range rolled down hill to the Exo Pro is comfy and extremely breathable.

Available on:

- ✓ Krypton R KRR030317
- ✓ Krypton F KRF030317
- ✓ Helium HUM030317
- ✓ Xenon² FF XFF030317
- ✓ Xenon² H XFF030317
- ✓ Xenon² SA XES030317
- ✓ Argon² AGN030317
- ✓ Neon² NSW030317
- ✓ Life F LFE030317
- ✓ Life R LRF030317




SIDEGUARDS AND ARMRESTS

Slim style carbon side guards

The slim style sideguards offer limited support and protection but are very lightweight for users who are not doing so much outside and have good trunk stability.

Available on:

- ✓ Krypton R KRR040101
- ✓ Krypton F KRF040101
- ✓ Helium HUM040101
- ✓ Xenon² FF XFF040101
- ✓ Xenon² H XFF040101
- ✓ Xenon² SA XES040101
- ✓ Argon² AGN040101



Straight side guard aluminium (without fender)

This sideguard offers a lot of support and protection from water and dirt from the tyres.

Available on:

- ✓ Krypton R KRR040102
- ✓ Krypton F KRF040102
- ✓ Helium HUM040102
- ✓ Xenon² FF XFF040102
- ✓ Xenon² H XFF040102
- ✓ Xenon² SA XES040102
- ✓ Argon² AGN040102
- ✓ Neon² NSW040102
- ✓ Life F LFE040102
- ✓ Life R LRF040102



Fender sideguard - aluminium

This sideguard offers the most support and protection against water and dirt from the wheels. However, it may not be suitable for users who grip the wheel in propulsion.

Available on:

- ✓ Krypton R KRR040103
- ✓ Krypton F KRF040103
- ✓ Helium HUM040103
- ✓ Xenon² FF XFF040103



Fender sideguard - aluminium - continued

✓ Xenon ² H	XFF040103
✓ Xenon ² SA	XES040103
✓ Argon ²	AGN040103
✓ Neon ²	NSW040103
✓ Life F	LFE040103
✓ Life R	LRF040103

Fender sideguard - carbon

Same as the fender sideguard (previous page) but in lightweight carbon fibre.

Available on:

✓ Krypton R	KRR040104
✓ Krypton F	KRF040104
✓ Helium	HUM040104
✓ Xenon ² FF	XFF040104
✓ Xenon ² H	XFF040104
✓ Xenon ² SA	XES040104
✓ Argon ²	AGN040104
✓ Neon ²	NSW040104



Plastic detachable sideguard (Quickie style)

These are a good compromise for someone who would like the protection of a sideguard but needs to remove them for transfers.

Available on:

✓ Krypton R	KRR040107
✓ Krypton F	KRF040107
✓ Helium	HUM040107
✓ Xenon ² FF	XFF040107
✓ Xenon ² H	XFF040107
✓ Xenon ² SA	XES040107
✓ Argon ²	AGN040107
✓ Neon ²	NSW040107
✓ Life F	LFE040107
✓ Life R	LRF040107



SIDEGUARDS AND ARMRESTS

Single post armrest - short armpad 25cm/33cm/36cm

This armrest offers support and protection with the ability to change the height of the armpad. The armpad is available in 25cm, 33cm and 36cm long.

25cm available on:

✓ Helium	HUM040050
✓ Xenon ² FF	XFF040050
✓ Xenon ² H	XFF040050
✓ Xenon ² SA	XES040050
✓ Argon ²	AGN040050
✓ Neon ²	NSW040050
✓ Life F	LFE040050
✓ Life R	LRF040050

33cm Available on:

✓ Helium	HUM040051
✓ Xenon ² FF	XFF040052
✓ Xenon ² H	XFF040052
✓ Xenon ² SA	XES040052
✓ Neon ²	NSW040052
✓ Life R	LRF040052

36cm Available on:

✓ Xenon ² FF	XFF040051
✓ Xenon ² H	XFF040051
✓ Xenon ² SA	XES040051
✓ Argon ²	AGN040051
✓ Neon ²	NSW040051
✓ Life F	LFE040051
✓ Life R	LRF040051




Single post, tool height adjustable armrest, depth adjustable (33cm)

Our lightest armrest offers the most support and protection with a built-in sideguard. It can be adjusted and moved forward/backward, as well as being height adjustable.

Available on:

✓ Krypton R	KRR040052
✓ Krypton F	KRF040052
✓ Helium	HUM040052
✓ Argon ²	AGN040052
✓ Life F	LFE040052
✓ Life R	LRF040052





Padded tubular armrest (swing- away)

This option is perfect for those who do not require an armrest all the time or just like to prop themselves up. They do not get in the way while propelling and can easily be swung-away or removed.

Available on:

- ✓ Krypton R KRR040117
- ✓ Krypton F KRF040117
- ✓ Helium HUM040108
- ✓ Xenon² FF XFF040117
- ✓ Xenon² H XFF040117
- ✓ Xenon² SA XES040117
- ✓ Argon² AGN040108
- ✓ Neon² NSW040117




Desk armrest with short armpad fixed (25cm)

This armrest offers a lot of protection, with a choice in pad lengths. The height of this armpad cannot be altered.

Available on:

- ✓ Xenon² FF XFF040001
- ✓ Xenon² H XFF040001
- ✓ Xenon² SA XES040001
- ✓ Neon² NSW040001
- ✓ Life F LFE040001



Similar options

Desk armrest with long armpad fixed (31cm)

Available on:

- ✓ Xenon² FF XFF040003
- ✓ Xenon² H XFF040003
- ✓ Xenon² SA XES040003
- ✓ Neon² NSW040003
- ✓ Life F LFE040003

Desk armrest without armpad

As above, without pad

Available on:

- ✓ Xenon² H XFF040005
- ✓ Xenon² SA XE040005
- ✓ Neon² NSW040005
- ✓ Life F LFE040005

Continues on following page...

SIDEGUARDS AND ARMRESTS

Similar options continued...

Desk armrest with short armpad height adj. (25cm)

As previous, but adjustable in height.

Available on:

- | | |
|-------------------------|-----------|
| ✓ Xenon ² FF | XFF040002 |
| ✓ Xenon ² H | XFF040002 |
| ✓ Xenon ² SA | XES040002 |
| ✓ Neon ² | NSW040002 |
| ✓ Life F | LFE040002 |

Desk armrest with long armpad height adjustable (31cm)

Available on:

- | | |
|-------------------------|-----------|
| ✓ Xenon ² FF | XFF040004 |
| ✓ Xenon ² H | XFF040004 |
| ✓ Xenon ² SA | XES040004 |
| ✓ Neon ² | NSW040004 |
| ✓ Life F | LFE040004 |

LEGRESTS AND FOOTRESTS

Lightweight tubular footrest

This is a simple lightweight solution. Good for users who do not require a lot of foot support.

Available on:

- ✓ Krypton R KRR050029
- ✓ Helium HUM050029
- ✓ Argon² AGN050029



Tubular footrest with flat cover

A small cover is attached to the tubular footrest to offer a little more support.

Available on:

- ✓ Krypton R KRR050030
- ✓ Helium HUM050030
- ✓ Argon² AGN050030



Performance footplate angle adjustable

This is our lightest footplate and offers good support and angle adjustment.

Available on:

- ✓ Krypton R KRR050125
- ✓ Krypton F KRF050059
- ✓ Helium HUM050125
- ✓ Xenon² FF XFF050059
- ✓ Xenon² H XFF050059
- ✓ Xenon² SA XES050059
- ✓ Argon² AGN050125
- ✓ Neon² NSW050059



LEGRESTS AND FOOTRESTS

Platform footplate lightweight angle and depth adjustable - carbon

This footplate can be adjusted forward and backward and angled. This offers maximum support and is available in composite, aluminium or carbon.

Available on:

- ✓ Krypton R KRR050026
- ✓ Helium HUM050026
- ✓ Argon² AGN050026



Similar options

Platform footplate lightweight angle and depth adjustable - aluminium

As previous, but in aluminium.

Available on:

- ✓ Krypton R KRR050024
- ✓ Helium HUM050024
- ✓ Argon² AGN050024
- ✓ Life R LRF050024

Platform footplates lightweight angle and depth adjustable composite

As previous, but in composite material.

Available on:

- ✓ Krypton R KRR050025
- ✓ Helium HUM050025
- ✓ Argon² AGN050025
- ✓ Life R LRF050025



Platform footplate lightweight angle & depth adjustable aluminium, high mounted

The high mount footplate flips back during transfer and is perfect for users with shorter legs.

Available on:

- ✓ Helium HUM050034
- ✓ Argon² AGN050034



Similar options

High mount platform footplates - composite footboard

As previous, but with a composite footboard.

Available on:

- ✓ Xenon² FF XFF050061
- ✓ Xenon² H XFF050061
- ✓ Xenon² SA XES050051
- ✓ Neon² NSW050051

Foot positioning plate

Foot positioning plates can be used to keep the feet in place and stop the feet sliding back off the footplate.

Available on:

- ✓ Krypton R KRR050127
- ✓ Krypton F KRF050127
- ✓ Helium HUM050127
- ✓ Xenon² FF XFF050127
- ✓ Xenon² H XFF050127
- ✓ Xenon² SA XES050127
- ✓ Argon² AGN050127
- ✓ Neon² NSW050127
- ✓ Life R LRF050127



LEGRESTS AND FOOTRESTS

Platform footrest, angle adjustable aluminium board

This flip-up footplate can be adjusted forward and backward and angled. This offers maximum support and is available in composite, aluminium or carbon.

Available on:

- ✓ Helium HUM050122
- ✓ Xenon² FF XFF050023
- ✓ Xenon² H XFF050023
- ✓ Xenon² SA XES050023
- ✓ Neon² NSW050023
- ✓ Life F LFE050102
- ✓ Life R LRF050102



Similar options

Platform footrest, angle adjustable composite board

This flip-up footplate can be adjusted forward and backward and angle. This offers maximum support and is available in composite, aluminium or carbon.

Available on:

- ✓ Krypton F KRF050131
- ✓ Xenon² FF XFF050131
- ✓ Xenon² H XFF050131
- ✓ Xenon² SA XES050131
- ✓ Neon² NSW050131
- ✓ Life F LFE050131
- ✓ Life R LRF050131

Platform footrest, angle adjustable carbon board

This flip-up footplate can be adjusted forward and backward and angle. This offers maximum support and is available in composite, aluminium or carbon.

Available on:

- ✓ Krypton F KRF050132
- ✓ Xenon² FF XFF050132
- ✓ Xenon² H XEFF050132
- ✓ Xenon² SA XES050132
- ✓ Neon² NSW050132



Auto-folding footrest angle and depth adjustable - carbon board

This auto flip-up footplate can be adjusted forward and backward and angled. This offers maximum support and is available in composite, aluminium or carbon. Automatically flips-up when the chair is folded.

Available on:

- ✓ Krypton F KRF050027
- ✓ Xenon² FF XFF050027
- ✓ Xenon² H XFF050027



Similar options

Auto-folding footrest angle & depth adjustable - composite board

As above, but in composite material.

Available on:


- ✓ Krypton F KRF050028
- ✓ Xenon² FF XFF050028
- ✓ Xenon² H XFF050028

Swing-away legrests - 70°

The Xenon legrests have a release catch that can be pushed or pulled horizontally for ease of use. This means that the trigger can be pushed or pulled in the direction that the legrest needs to swing in or out. The Xenon is one of our high-end manual chairs. The Xenon hangers are our stiffest most robust offering and are available across our other chairs.

Available on:

- ✓ Xenon² SA XES050002
- ✓ Neon² NSW050002
- ✓ Life F LFE050002
- ✓ Life R LRF050002



LEGRESTS AND FOOTRESTS

Similar options

Swing-away legrests - 80°

The Xenon legrests have a release catch that can be pushed or pulled horizontally for ease of use. This means that the trigger can be pushed or pulled in the direction that the legrest needs to swing in or out. The Xenon is one of our high end manual chairs. The Xenon hangers are our stiffest most robust offering and are available across our other chairs.

Available on:

- ✓ Xenon² SA XES050001
- ✓ Neon² NSW050001
- ✓ Life F LFE050001
- ✓ Life R LRF050001

Helix swing-away legrests - 70°

The Helix hangers have more composite built into the design so are not as stiff as the Xenon equivalents. The trigger is pulled toward the user whilst in the chair and may be better for users that have limited hand function.

Available on:

- ✓ Life F LFE050002
- ✓ Life R LRF050002



Similar options

HeliX swing-away legrests - 80°

As above but with an 80° angle.

Available on:

- ✓ Life F LFE050001
- ✓ Life R LRF050001



HeliX type Elevating legrests - pair

Elevated legrests (ELRs) can be used either as a single or as a pair when there is a requirement to raise the users legs.

Available on:

- ✓ Life F Pick both

Elevating legrests - left

Available on:

- ✓ Life F LFE050003

Elevating legrests - right

Available on:

- ✓ Life F LFE050004

Elevating legrests - pair

Elevated legrests (ELRs) can be used either as a single or as a pair when there is a requirement to raise the users legs.

Available on:

- ✓ Xenon² SA XES050011
- ✓ Neon² NSW050011
- ✓ Life R LRF050011



Similar options

Elevating legrests - left

Available on:

- ✓ Xenon² SA XES050009
- ✓ Neon² NSW050009
- ✓ Life R LRF050009

Elevating legrests - right

Available on:

- ✓ Xenon² SA XES050010
- ✓ Neon² NSW050010
- ✓ Life R LRF050010

LEGRESTS AND FOOTRESTS

Angle adjustable footplates (composite)

These footplates can be flipped-up to create space at the front of the chair to transfer and have angle adjustment to accommodate the user. The depth can be adjusted slightly by swapping the footplates over (right to left).

Available on:

- ✓ Xenon² SA XES050020
- ✓ Argon² AGN050020
- ✓ Neon² NSW050020
- ✓ Life F LFE050020
- ✓ Life R LRF050020

Divided footplates, aluminium

These footplates are more hard-wearing than the composite equivalents and can also be angle and depth adjusted.

Available on:

- ✓ Xenon² SA XES050022
- ✓ Argon² AGN050022
- ✓ Neon² NSW050022
- ✓ Life F LFE050050
- ✓ Life R LRF050050



Amputee pad - right

This pad supports the amputated leg, this can be ordered as left/right or pair (see opposite). It is also worth considering a passive or amputee axle to avoid the chair being too tippy.

Available on:

- ✓ Xenon² SA XES050013
- ✓ Life F LFE050013





Similar options

Amputee pad - left

Available on:

- ✓ Xenon² SA XES050014
- ✓ Life F LFE050014

Amputee pad - pair

Available on:

- ✓ Xenon² SA XES050012
- ✓ Life F LFE050012


FRAME VERSIONS

Frame

Frame angle

Frame angle is taken from the horizontal (grey line) to the front frame, not the inside of the frame.

This means that a 90° frame will have a smaller overall footprint and will cater for those with closed knee angles; whereas a 75° frame would make the overall chair a little longer and would cater for those with more open knee angles.



Frame

Long frame

The frame is a couple of centimetres longer than the seat canvas as per the traditional Helium.

Available on:

- ✓ Krypton R KRR010207
- ✓ Helium HUM010207



Similar options

Standard frame

Available on:

- ✓ Krypton R KRR010206
- ✓ Helium HUM010206

Similar options (continued from previous page)

Frame angle 100° (UK 80°)

Available on:

- ✓ Krypton R KRR010012
- ✓ Krypton F KRF010012
- ✓ Helium HUM010012
- ✓ Xenon² FF XES010012
- ✓ Xenon² H XES010020
- ✓ Argon² AGN010012

Frame angle 92° (UK 88°)

Available on:

- ✓ Krypton R KRR010013
- ✓ Krypton F KRF010013
- ✓ Helium HUM010013
- ✓ Xenon² FF XES010013
- ✓ Xenon² H XES010021
- ✓ Argon² AGN010013

Frame angle 85°

Available on:

- ✓ Life F LFE010011
- ✓ Life R LRF010011

Frame angle 75°

Available on:

- ✓ Life F LFE010010
- ✓ Life R LRF010010

Frame angle 85° (2cm inset)

Available on:

- ✓ Life F LFE010016

Frame angle 85° (-2cm inset)

Available on:

- ✓ Life R LRF010016


FRAME VERSIONS

Inset: Standard (0cm)

This frame does not taper at the front and provides the maximum space available at the front of the chair for legs and feet.

Available on:

- ✓ Helium HUM010100
- ✓ Argon² AGN010100



Similar options

Inset: 1cm

This frame tapers slightly at the front of the chair by 1cm each side. This offers a sportier, sleeker look without encroaching on the space of the footplate and legs.

Available on:

- ✓ Krypton R KRR010101
- ✓ Krypton F KRF010012
- ✓ Helium HUM010101
- ✓ Xenon² FF XFF010101
- ✓ Xenon² H XFF010101

Inset: 2cm

This frame has a more aggressive taper, 2cm each side that will start to reduce the front footprint of the chair.

Available on:

- ✓ Life R LRF010016
- ✓ Life F LFE010016




Inset: (3cm) only for frame angle 92°

This frame has our maximum taper of 3cm each side on the Helium which offers a very sleek look and small front footprint. Please bear in mind that the footplate will be 12cm shorter than the width of the chair, always worth checking what footwear the user likes to wear (especially if they wear bulky footwear) particularly on narrow frame widths.

Available on:

- ✓ Krypton R KRR010103
- ✓ Helium HUM010103
- ✓ Argon² AGN010103




Castor house position wide

This is the traditional setting for the Helium. The castor wheels are very aggressively set outwards at the front of the chair and provide a lot of forward stability.

Available on:

- ✓ Krypton R KRR010301
- ✓ Helium HUM010301



Similar options

Castor house position narrow

The narrow position brings the front castor wheels in by 3cm. This decreases the front frame footprint and brings the castor wheels inside the rear wheels. This can make the front of the chair more tippy when leaning forwards or transferring. To ensure the front stability of the chair we recommend this setting in combination with our 'long frame' option. This is also our most energy efficient set-up.

Available on:

- ✓ Krypton R KRR010302
- ✓ Helium HUM010302

CENTRE OF GRAVITY

Centre of gravity


The centre of gravity measurement is taken from the front of the backposts forward or backward.

In the diagram;

Position A would be considered an active setting, for example 10cm.

Position B may be classed as a passive or amputee setting, in this case this may -1cm.

Position C gives an example of a handbike setting, very far behind the backposts, the measurements for these are not generally given but are usually around -6cm.



Helium Standard centre of gravity

The axle position is in front of the backposts and can be adjusted forward or backward. The amount of adjustment is dependant on other options that are chosen and the length of the frame. Please note that there is only +/-1cm of COG adjustment on the Helium Pro.

Simply select this option on the order form.



Helium passive centre of gravity (amputee)

This axle position is in line and behind the backposts to ensure the chair is not too tippy and is a good feature for those that are being pushed the majority of the time, or for amputees.

Simply select this option on the order form.

Helium Handbike

The handbike axle places the wheels behind the backposts so that the handbike wheel gains more traction when attached. The handbike axle for the Helium and Life R are separate to the everyday axle and can be kept on the chair during everyday use.

Available on:

- ✓ Helium HUM010500



Adapter Lightweight

The lightweight axle is set in front of the backposts but only has one position which can be adjusted but takes more time. This is the lightest option available.

Available on:

- ✓ Krypton F KRF070701
- ✓ Xenon² FF XFF070701
- ✓ Xenon² H XFF070701
- ✓ Xenon² SA XES070701
- ✓ Neon² NSW070701



CENTRE OF GRAVITY

Adapter Standard

This axle position is inline and behind the backposts to ensure the chair is not too tippy and is a good feature for those that are being pushed the majority of the time or for amputees.

Available on:

- ✓ Krypton F KRF070703
- ✓ Xenon² FF XFF070703
- ✓ Xenon² H XFF070703
- ✓ Xenon² SA XES070703
- ✓ Neon² NSW070703



Argon - Standard

The axle position is in front of the backposts and can be adjusted forward or backward. The amount of adjustment is dependant on other options that are chosen and the length of the frame.

Simply select this option on the order form.



Argon Passive/Amputee

This axle position is in-line and behind the backposts to ensure the chair is not too tippy and is a good feature for those that are being pushed the majority of the time, or for amputees.

Available on:

- ✓ Argon² AGN070500



Argon Handbike Axle

The handbike axle places the wheels behind the backposts so that the handbike wheel gains more traction when attached. The handbike axle for the Helium and Life R are separate to the everyday axle and can be kept on the chair during everyday use.

Available on:

✓ Argon² AGN010500



Life/Life R - Standard

The Life and Life R standard axle is in front of the backposts and all of the COG options available on the order form are available for adjustment later on.

Simply select this option on the order form.



Life R - Handbike Axle

The handbike axle places the wheels behind the backposts so that the handbike wheel gains more traction when attached. The handbike axle for the Helium and Life R are separate to the everyday axle and can be kept on the chair during everyday use.

Available on:

✓ Life R LRF010500



FRAME ACCESSORIES

Suspension up to 70kg

Suspension is built in to the rear axle and can dampen the ride across rough surfaces. This may be a good option to reduce spasm, however, this does add weight to the chair.

Available on:

- ✓ Krypton R KRR090031
- ✓ Helium HUM090031



Similar options

Suspension up to 120kg

Suspension is built into the rear axle and can dampen the ride across rough surfaces. This maybe a good option to reduce spasm, however, this does add weight to the chair.

Available on:

- ✓ Krypton R KRR090032
- ✓ Helium HUM090032

Tip assist (step tube) Left

Tip assist tubes are useful for carers to step on when tipping the chair back when going over kerbs and thresholds.

Available on:

- ✓ Krypton R KRR090001
- ✓ Krypton F KRF090001
- ✓ Xenon² FF XFF090001
- ✓ Xenon² H XFF090001
- ✓ Xenon² SA XES090001
- ✓ Neon² NSW090001
- ✓ Argon² AGN090001
- ✓ Life F LFE090001
- ✓ Life R LRF090001



Similar options

Tip assist (step tube) Right

Available on:

- | | |
|-------------------------|-----------|
| ✓ Krypton R | KRR090002 |
| ✓ Krypton F | KRF090002 |
| ✓ Helium | HUM090002 |
| ✓ Xenon ² FF | XFF090002 |
| ✓ Xenon ² H | XFF090002 |
| ✓ Xenon ² SA | XES090002 |
| ✓ Argon ² | AGN090002 |
| ✓ Neon ² | NSW090002 |
| ✓ Life F | LFE090002 |
| ✓ Life R | LRF090002 |

Tip assist (step tube) Pair

Available on:

- | | |
|-------------------------|-----------|
| ✓ Krypton R | KRR090003 |
| ✓ Krypton F | KRF090003 |
| ✓ Helium | HUM090003 |
| ✓ Xenon ² FF | XFF090003 |
| ✓ Xenon ² H | XFF090003 |
| ✓ Xenon ² SA | XES090003 |
| ✓ Argon ² | AGN090003 |
| ✓ Neon ² | NSW090003 |
| ✓ Life F | LRE090003 |
| ✓ Life R | LRF090003 |

Crutch holder

The crutch holder is mounted at the rear of the frame so crutches can be stowed while in the chair.
(Requires a tip assist to mount).

Available on:

- | | |
|-------------------------|-----------|
| ✓ Krypton R | KRR090009 |
| ✓ Krypton F | KRF090009 |
| ✓ Helium | HUM090009 |
| ✓ Xenon ² FF | XFF090009 |
| ✓ Xenon ² H | XFF090009 |
| ✓ Xenon ² SA | XES090009 |
| ✓ Argon ² | AGN090009 |
| ✓ Neon ² | NSW090009 |
| ✓ Life F | LFE090009 |
| ✓ Life R | LRF090009 |



FRAME ACCESSORIES

Anti-tip tubes swing-away - left

These anti-tip tubes can be swung inward under the chair when not required or when being pushed.

Available on:

✓ Krypton R	KRR090004
✓ Krypton F	KRF090004
✓ Helium	HUM090004
✓ Xenon ² FF	XFF090004
✓ Xenon ² H	XFF090004
✓ Xenon ² SA	XES090004
✓ Argon ²	AGN090004
✓ Neon ²	NSW090004
✓ Life F	LFE090004
✓ Life R	LRF090004



Similar options

Anti-tip tubes swing-away - right

Available on:

✓ Krypton R	KRR090005
✓ Krypton F	KRF090005
✓ Helium	HUM090005
✓ Xenon ² FF	XFF090005
✓ Xenon ² H	XFF090005
✓ Xenon ² SA	XES090005
✓ Argon ²	AGN090005
✓ Neon ²	NSW090005
✓ Life F	LFE090005
✓ Life R	LRF090005

Anti-tip tubes swing-away - pair

Available on:

✓ Krypton R	KRR090006
✓ Krypton F	KRF0900056
✓ Helium	HUM090006
✓ Xenon ² FF	XFF090006
✓ Xenon ² H	XFF090006
✓ Xenon ² SA	XES090006
✓ Argon ²	AGN090006
✓ Neon ²	NSW090006
✓ Life F	LFE090006
✓ Life R	LRF090006

Quickie style anti-tips - pair

These anti-tips can be detached or rotated upwards when not required. These are our lightest option.

Available on:

✓ Krypton R	KRR090007
✓ Krypton F	KRF0900050
✓ Helium	HUM090007
✓ Xenon ² FF	XFF0900050
✓ Xenon ² H	XFF0900050
✓ Xenon ² SA	XES0900050





Continued from previous page

- ✓ Argon² AGN090007
- ✓ Life F LIF090036
- ✓ Life R LIR090036

Anti-tip flip-up (user operated) - Left

This anti-tip is a great compromise for a user that does not require an anti-tip all the time. The user must have the trunk support to reach behind and under the chair whilst seated.

Available on:

- ✓ Krypton R KRR090033
- ✓ Helium HUM090008
- ✓ Argon² AGN090033



Similar options

Anti-tip flip-up (user operated) - right

Available on:

- ✓ Krypton R HUM090008
- ✓ Helium HUM090008
- ✓ Argon² AGN090008

Sport anti-tip, central mount, detachable

This anti-tip is perfect for those who would like to give non-contact wheelchair sports a try (like badminton or dancing). It is mounted at the centre of the axle and can be adjusted to the correct distance away from the floor, it can also be detached when not in use.

Available on:

- ✓ Helium HUM090035
- ✓ Argon² AGN090035



FRAME ACCESSORIES

Transit wheels

Transit wheels attach at the rear of the frame and are commonly used on planes and overseas through tight spaces and doorways.

Available on:

- | | |
|-------------------------|-----------|
| ✓ Helium | HUM090010 |
| ✓ Xenon ² FF | XFF090010 |
| ✓ Xenon ² H | XFF090010 |
| ✓ Xenon ² SA | XES090010 |
| ✓ Argon ² | AGN090010 |
| ✓ Neon ² | NSW090010 |
| ✓ Life F | LFE090010 |
| ✓ Life R | LRF090010 |



Frame protection

These covers can protect the framework at the front of the chair from everyday scrapes and scratches.

Available on:

- | | |
|-------------------------|-----------|
| ✓ Krypton R | KRR090016 |
| ✓ Krypton F | KRF090016 |
| ✓ Helium | HUM090016 |
| ✓ Xenon ² FF | XFF090016 |
| ✓ Xenon ² H | XFF090016 |
| ✓ Xenon ² SA | XES090016 |
| ✓ Argon ² | AGN090016 |
| ✓ Neon ² | NSW090016 |




FORKS

Carbotecture fork

These are our lightest forks and are made from a hybrid of composite and carbon fibre. They are very strong and offer a little flex that can dampen the ride slightly.

Available on:

- | | |
|-------------------------|-----------|
| ✓ Krypton R | KRR080013 |
| ✓ Krypton F | KRF080013 |
| ✓ Helium | HUM080007 |
| ✓ Xenon ² FF | XFF080013 |
| ✓ Xenon ² H | XFF080013 |
| ✓ Xenon ² SA | XES080013 |
| ✓ Argon ² | AGN080013 |
| ✓ Neon ² | NSW080013 |



Single arm fork

Single arm forks decrease the front frame footprint slightly and add an element of style to the front of the chair.

Available on:

- | | |
|-------------------------|-----------|
| ✓ Krypton R | KRR080008 |
| ✓ Krypton F | KRF080008 |
| ✓ Helium | HUM080008 |
| ✓ Xenon ² FF | XFF080008 |
| ✓ Xenon ² H | XFF080008 |
| ✓ Xenon ² SA | XEOS80008 |
| ✓ Neon ² | NSW080008 |



Frog Legs fork

Frog legs have built in suspension and are great for decreasing spasm, but do add weight.

Available on:

- | | |
|-------------------------|-----------|
| ✓ Krypton R | KRR080007 |
| ✓ Krypton F | KRF080007 |
| ✓ Helium | HUM080007 |
| ✓ Xenon ² FF | XFF080007 |
| ✓ Xenon ² H | XFF080007 |
| ✓ Xenon ² SA | XEOS80007 |
| ✓ Argon ² | AGN080007 |



FORKS

Castor fork for power add-on

Made out of aluminium to give strength to the chair when used with power add-ons.

Available on:

✓ Helium	HUM090093
✓ Xenon ² FF	XFF090093
✓ Xenon ² H	XFF090093
✓ Xenon ² SA	XES090093
✓ Argon ²	AGN090093
✓ Neon ²	NSW090093



CASTOR WHEELS

Solid

Solid castor wheels will last longer than most but can make the ride feel harder over rough surfaces.

Available on:

- ✓ Krypton R KRR080301
- ✓ Krypton F KRF080301
- ✓ Helium HUM080301
- ✓ Xenon² FF XFF080301
- ✓ Xenon² H XFF080301
- ✓ Xenon² SA XES080301
- ✓ Argon² AGN080301
- ✓ Neon² NSW080301
- ✓ Life F LFE080301
- ✓ Life R LRF080301



Soft

Soft castors may wear faster but do provide a smoother ride over rough surfaces. Also feature a wider profile.

Available on:

- ✓ Krypton R KRR080302
- ✓ Krypton F KRF080302
- ✓ Helium HUM080302
- ✓ Xenon² FF XFF080302
- ✓ Xenon² H XFF080302
- ✓ Xenon² SA XES080302
- ✓ Argon² AGN080302
- ✓ Neon² NSW080302
- ✓ Life F LFE080302
- ✓ Life R LRF080302



Pneumatic

Pneumatic castors can provide a smooth ride across rough surfaces and will climb over small obstacles much better than others although they do need to be maintained and inflated.

Available on:

- ✓ Xenon² SA XES080305
- ✓ Neon² NSW080305
- ✓ Life F LFE080305
- ✓ Life R LRF080305



CASTOR WHEELS

Solid light up

The solid light ups will last longer and ensure the chair looks eye catching!

Available on:

- ✓ Krypton R KRR080303
- ✓ Krypton F KRF080303
- ✓ Helium HUM080303
- ✓ Xenon² FF XFF080303
- ✓ Xenon² H XFF080303
- ✓ Xenon² SA XES080303
- ✓ Argon² AGN080303
- ✓ Neon² NSW080303



Solid with aluminium rim

These castor wheels look sleek and will last longer but do provide a harder ride over rough surfaces.

Available on:

- ✓ Krypton R KRR080304
- ✓ Krypton F KRF080304
- ✓ Helium HUM080304
- ✓ Xenon² FF XFF080304
- ✓ Xenon² H XFF080304
- ✓ Xenon² SA XES080304
- ✓ Argon² AGN080304
- ✓ Neon² NSW080304



Soft with aluminium rim

These castor wheels look sleek and will give a softer ride over rougher surfaces but may wear faster than harder castors. Also feature a wider profile.

Available on:

- ✓ Krypton R KRR080306
- ✓ Krypton F KRF080306
- ✓ Helium HUM080306
- ✓ Xenon² FF XFF080306
- ✓ Xenon² H XFF080306
- ✓ Xenon² SA XES080306
- ✓ Argon² AGN080306
- ✓ Neon² NSW080306



REAR WHEELS

Universal Cross spoked

Standard wheels that require minimal maintenance.

Available on:

- ✓ Helium HUM070001
- ✓ Xenon² FF XFF070001
- ✓ Xenon² H XFF070001
- ✓ Xenon² SA XES070001
- ✓ Argon² AGN070001
- ✓ Neon² NSW070001
- ✓ Life F LFE070001
- ✓ Life R LRF070001



Design wheel

Stylish high polished aluminium.

Available on:

- ✓ Helium HUM070002
- ✓ Xenon² FF XFF070002
- ✓ Xenon² H XFF070002
- ✓ Xenon² SA XES070002
- ✓ Argon² AGN070002
- ✓ Neon² NSW070002
- ✓ Life F LFE070002
- ✓ Life R LRF070002



Lightweight wheel (silver spokes)

Provide 1.5kg weight saving over a pair of straight spoke wheels.

Not available with solid tyres.

Available on:

- ✓ Krypton R KRR070003
- ✓ Krypton F KRF070003
- ✓ Helium HUM070003
- ✓ Xenon² FF XFF070003
- ✓ Xenon² H XFF070003
- ✓ Xenon² SA XES070003
- ✓ Argon² AGN070003
- ✓ Neon² NSW070003
- ✓ Life F LFE070003
- ✓ Life R LRF070003



REAR WHEELS

Proton wheel (super lightweight) (black spokes)

These are our lightest wheels and have a sleek look with black wheels and hub and proton branding.

Available on:

✓ Krypton R	KRR070004
✓ Krypton F	KRF070004
✓ Helium	HUM070004
✓ Xenon ² FF	XFF070004
✓ Xenon ² H	XFF070004
✓ Xenon ² SA	XES070004
✓ Argon ²	AGN070004
✓ Neon ²	NSW070004



Spinergy wheels, 18 spokes

Made with double-walled custom design aluminium rim for strength. Very little maintenance required.

Available on:

✓ Krypton R	KRR070005	✓ Life F	LFE070005
✓ Krypton F	KRF070005	✓ Life R	LRF070005
✓ Helium	HUM070005		
✓ Xenon ² FF	XFF070005		
✓ Xenon ² H	XFF070005		
✓ Xenon ² SA	XES070005		
✓ Argon ²	AGN070005		
✓ Neon ²	NSW070005		



Spinergy wheels, 12 spokes

Made with double-walled custom design aluminium rim for strength. Very little maintenance required.

Available on:

✓ Krypton R	KRR070006 (silver) & KRR070007 (black)
✓ Krypton F	KRF070006 (silver) & KRF070007 (black)
✓ Helium	HUM070006 (silver) & HUM070007 (black)
✓ Xenon ² FF	XFF070006 (silver) & XFF070007 (black)
✓ Xenon ² H	XFF070006 (silver) & XFF070007 (black)
✓ Xenon ² SA	XES070006 (silver) & XES070007 (black)
✓ Neon ²	NSW070006 (silver) & NSW070007 (black)
✓ Argon ²	AGN070006 (silver)





Mountain bike wheel

This wheel comes complete with a mountain bike tyre ready for off road use.

Available on:

- | | |
|-------------------------|-----------|
| ✓ Helium | HUM070009 |
| ✓ Xenon ² FF | XFF070009 |
| ✓ Xenon ² H | XFF070009 |
| ✓ Xenon ² SA | XES070009 |
| ✓ Argon ² | AGN070009 |
| ✓ Neon ² | NSW070009 |



Similar options

Additional mountain bike wheel (pair)

Available on:

- | | |
|-------------------------|-----------|
| ✓ Helium | HUM070010 |
| ✓ Xenon ² FF | XFF070010 |
| ✓ Xenon ² H | XFF070010 |
| ✓ Xenon ² SA | XES070010 |
| ✓ Argon ² | AGN070010 |
| ✓ Neon ² | NSW070010 |

Right Run

This air filled tyre offers a cushioned ride and a slick tread for low rolling resistance. Available in 22", 24" (Schwalbe right run) and 26" (slick tyre).

Available on:

- | | |
|-------------------------|-----------|
| ✓ Krypton R | KRR070101 |
| ✓ Krypton F | KRF070101 |
| ✓ Helium | HUM070101 |
| ✓ Xenon ² FF | XFF070101 |
| ✓ Xenon ² H | XFF070101 |
| ✓ Xenon ² SA | XES070101 |
| ✓ Argon ² | AGN070101 |
| ✓ Neon ² | NSW070101 |
| ✓ Life F | LFE070101 |
| ✓ Life R | LRF070101 |



REAR WHEELS

Down Town Two Grip

At each push and particularly when braking, the hands constant contact with tyres can result in irritation and chafing. The Schwalbe Downtown with its patented 2 grip solution consists of completely smooth sidewalls, allowing hands to glide past the tyres without irritation.

Available on:

- | | |
|-------------------------|-----------|
| ✓ Krypton R | KRR070106 |
| ✓ Krypton F | KRF070106 |
| ✓ Helium | HUM070106 |
| ✓ Xenon ² FF | XFF070106 |
| ✓ Xenon ² H | XFF070106 |
| ✓ Xenon ² SA | XES070106 |
| ✓ Argon ² | AGN070106 |
| ✓ Neon ² | NSW070106 |



Solid Tyres

Solid tyres can provide a smooth ride over smooth surfaces, however, over rough surfaces can be very hard and it can also be difficult to climb kerbs. They are also maintenance free as they do not require inflation.

Available on:

- | | |
|-------------------------|-----------|
| ✓ Krypton R | KRR070102 |
| ✓ Krypton F | KRF070102 |
| ✓ Helium | HUM070102 |
| ✓ Xenon ² FF | XFF070102 |
| ✓ Xenon ² H | XFF070102 |
| ✓ Xenon ² SA | XES070102 |
| ✓ Argon ² | AGN070102 |
| ✓ Neon ² | NSW070102 |
| ✓ Life F | LFE070102 |
| ✓ Life R | LRF070102 |



Marathon Plus Evolution

Extremely effective at preventing defects. The insert, made from highly elastic special rubber, offers durable resistance to penetration punctures. Combines comfort and easy rolling of a pneumatic tyre with the puncture resistance of solid tyres.

Available on:

✓ Krypton R	KRR070104
✓ Krypton F	KRF070104
✓ Helium	HUM070104
✓ Xenon ² FF	XFF070104
✓ Xenon ² H	XFF070104
✓ Xenon ² SA	XES070104
✓ Argon ²	AGN070104
✓ Neon ²	NSW070104
✓ Life F	LIF080513
✓ Life R	LIR070104



Schwalbe One

Available on:

✓ Krypton R	KRR070107
✓ Krypton F	KRF070107
✓ Helium	HUM070107
✓ Xenon ² FF	XFF070104
✓ Xenon ² H	XFF070104
✓ Xenon ² SA	XSE070104
✓ Argon ²	AGN070107
✓ Neon ²	NSW070107
✓ Life F	LFE070107
✓ Life R	LRF070107



HANDRIMS

Supergrip

A rubberised handrim that offers more grip and a solid feel.

Available on:

✓ Krypton R	KRR070207
✓ Krypton F	KRF070207
✓ Helium	HUM070207
✓ Xenon ² FF	XFF070207
✓ Xenon ² H	XFF070207
✓ Xenon ² SA	XES070207
✓ Argon ²	AGN070207
✓ Neon ²	NSW070207
✓ Life F	LFE070207
✓ Life R	LRF070207



Max Grepp

This option fills the gap between the handrim and wheel with a soft rubber making it easier to propel for those with less hand function.

Available on:

✓ Krypton R	KRR070208
✓ Krypton F	KRF070208
✓ Helium	HUM070208
✓ Xenon ² FF	XFF070208
✓ Xenon ² H	XFF070208
✓ Xenon ² SA	XES070208
✓ Argon ²	AGN070208
✓ Neon ²	NSW070208
✓ Life F	LFE070208
✓ Life R	LRF070208



Handrim cover

A rubber cover that is fitted over the handrim that can also be removed.

Available on:

✓ Krypton R	KRR070250
✓ Krypton F	KRF070250
✓ Helium	HUM070250
✓ Xenon ² FF	XFF070250
✓ Helium	HUM070250
✓ Xenon ² FF	XFF070250
✓ Xenon ² H	XFF070250
✓ Xenon ² SA	XES070250
✓ Argon ²	AGN070250
✓ Neon ²	NSW070250
✓ Krypton F	LFE070250
✓ Krypton R	LRF070250



Surge LT

The lightweight handrim offers gripton strips which improve traction and the oval shape fits nicely to the shape of the hand

Available on:

✓ Krypton R	KRR070211
✓ Krypton F	KRK070211
✓ Helium	HUM070211
✓ Xenon ² FF	XFF070211
✓ Xenon ² H	XFF070211
✓ Xenon ² SA	XES070211
✓ Argon ²	AGN070211
✓ Neon ²	NSW070211
✓ Life F	LFE070211
✓ Life R	LRF070211



HANDRIMS

Handrim mounting - Close 17mm

This is the standard setting that most users will be used to. It is now possible to have the handrim mounted closer at 7mm and spaced out at 27mm.

Available on:

✓ Krypton R	KRR070317
✓ Krypton F	KRF070317
✓ Helium	HUM070317
✓ Xenon ² FF	XFF070317
✓ Xenon ² H	XFF070317
✓ Xenon ² SA	XES070317
✓ Argon ²	AGN070317
✓ Neon ²	NSW070317
✓ Life F	LFE070317
✓ Life R	LRF070317



Similar options

Handrim mounting - Extra close 7mm

Available on:

✓ Krypton R	KRR070307
✓ Krypton F	KRF070307
✓ Helium	HUM070307
✓ Xenon ² FF	XFF070307
✓ Xenon ² H	XFF070307
✓ Xenon ² SA	XES070307
✓ Argon ²	AGN070307
✓ Neon ²	NSW070307
✓ Life F	LFE070307
✓ Life R	LRF070307

Handrim mounting - Wide 27mm

Available on:

✓ Krypton R	KRR070327
✓ Krypton F	KRF070327
✓ Helium	HUM070327
✓ Xenon ² FF	XFF070327
✓ Xenon ² H	XFF070327
✓ Xenon ² SA	XES070327
✓ Argon ²	AGN070327
✓ Neon ²	NSW070327
✓ Life F	LFE070327
✓ Life R	LRF070327

SPOKEGUARDS



Princess



Dancing Bear



Jungle Fever



Happy Lion



Water World



Randolf

Spokeguards are available across our full range of manual chairs, just select them on the order form.

SPOKEGUARDS



Monkeys



Pirates



Feather



Oriental



Racing



Dirt Race

Spokeguards are available across our full range of manual chairs, just select them on the order form.



Suphero



Graffiti



Chrome Wheel



Hot Wheels



Basketball




Basketball 2

Spokeguards are available across our full range of manual chairs, just select them on the order form.


SPOKEGUARDS



Pool



Flying Water



Space



Silk



Abstract



Creative

Spokeguards are available across our full range of manual chairs, just select them on the order form.



Tiger



Zebra



Cheese



Bibbo



Galactica



Rotation

Spokeguards are available across our full range of manual chairs, just select them on the order form.

SPOKEGUARDS



Fantasy



Graphic

Spokeguards are available across our full range of manual chairs, just select them on the order form.

WHEEL LOCKS

Standard wheel lock (Active wheel lock)

Spring loaded for ease of use, when engaged the lever lies out of the way making side transfers easy.

Available on:

✓ Krypton R	KRR060001
✓ Krypton F	KRF060001
✓ Helium	HUM060001
✓ Xenon ² FF	XFF060001
✓ Xenon ² H	XFF060001
✓ Xenon ² SA	XES060001
✓ Argon ²	AGN060001
✓ Neon ²	NSW060001
✓ Life F	LFE060001
✓ Life R	LRF060001




Knee lever wheel lock

Spring-loaded for ease of use and mounted at knee level for easy access, this wheel lock is suitable for most users.

Available on:

✓ Krypton R	KRR060002
✓ Krypton F	KRF060002
✓ Helium	HUM060002
✓ Xenon ² FF	XFF060002
✓ Xenon ² H	XFF060002
✓ Xenon ² SA	XES060002
✓ Argon ²	AGN060002
✓ Neon ²	NSW060002
✓ Life F	LFE060002
✓ Life R	LRF060002



WHEEL LOCKS

Extension for knee lever wheel lock, slip-on, swing-away

This extension handle makes operating the knee lever wheel lock easier for users with restricted movement.

Available on:

- | | |
|-------------------------|-----------|
| ✓ Krypton R | KRR060010 |
| ✓ Krypton F | KRF060010 |
| ✓ Helium | HUM060010 |
| ✓ Xenon ² FF | XFF060010 |
| ✓ Xenon ² H | XFF060010 |
| ✓ Xenon ² SA | XES060010 |
| ✓ Argon ² | AGN060010 |
| ✓ Neon ² | NSW060010 |
| ✓ Life F | LFE060010 |
| ✓ Life R | LRF060010 |




Compact wheel lock

This wheel lock is lightweight and is less likely to get in the way of propulsion. It may not be suitable for users that aggressively transfer out of the chair.

Available on:

- | | |
|-------------------------|-----------|
| ✓ Krypton R | KRR060003 |
| ✓ Krypton F | KRF060003 |
| ✓ Helium | HUM060003 |
| ✓ Xenon ² FF | XFF060003 |
| ✓ Xenon ² H | XFF060003 |
| ✓ Xenon ² SA | XES060003 |
| ✓ Argon ² | AGN060003 |
| ✓ Neon ² | NSW060003 |
| ✓ Life F | LFE060003 |
| ✓ Life R | LRF060003 |





Lightweight compact

This wheel lock is our lightest brake. It is also less likely to get in the way of propulsion however there is less adjustment for those who may wish to change their COG position and it may not be suitable for users that aggressively transfer out the chair.

Available on:

- ✓ Krypton R KRR060004
- ✓ Helium HUM060004



One arm wheel lock -left

Uses combined compact brakes on both sides of the chair with one side being operational. This brake is suitable for those with bias strength, one arm, or in some sports such as archery or shooting.

Available on:

- ✓ Krypton R KRR060006
- ✓ Helium HUM060006
- ✓ Argon² AGN060006
- ✓ Neon² NSW060015
- ✓ Life F LFE060015



Similar options

One arm wheel lock right

Available on:

- ✓ Krypton R KRR060007
- ✓ Helium HUM060007
- ✓ Argon² AGN060007
- ✓ Neon² NSW060016
- ✓ Life F LFE060016

ACCESSORIES

Tool kit

Contains the tools required for the adjustments available on the chair.

Available on:

✓ Krypton R	KRR090000	✓ Life F	LFE090000
✓ Krypton F	KRF090000	✓ Life R	LRF090000
✓ Helium	HUM090000		
✓ Xenon ² FF	XFF090000		
✓ Xenon ² H	XFF090000		
✓ Xenon ² SA	XES090000		
✓ Argon ²	AGN090000		
✓ Neon ²	NSW090000		



Snoll on

A fifth wheel attachment that fits to the footplate of the chair and raises the castor wheels off the floor. This is intended for use outdoor or prolonged use over carpet. The Snoll on has a single directional wheel and does not turn, if you require a wheel attachment that rotates then please contact our Built-4-Me department who will work with you to find a solution.

Available on:

✓ Krypton R	KRR090024
✓ Krypton F	KRF090024
✓ Helium	HUM070205
✓ Argon ²	AGN090017



High pressure pump

Most of our pneumatic tyres are high pressure cycling tyres and can be pressurised up to 140 psi (please refer to the recommended amount on the tyre). The high pressure pump can fill the tyre to this level although some users may require assistance to do so.

Available on:

✓ Krypton R	KRR090024	✓ Argon ²	AGN090024
✓ Krypton F	KRF090024	✓ Neon ²	NSW090024
✓ Helium	HUM090024	✓ Life F	LFE090024
✓ Xenon ² FF	XFF090024	✓ Life R	LRF090024
✓ Xenon ² H	XFF090024		
✓ Xenon ² SA	XES090024		



Air pump

A simple hand pump that can keep the tyre inflated to a reasonable pressure but is very difficult to reach the true pressure for the chair. Useful for the user to keep with them to keep punctures inflated while out and about until they can be repaired.

Available on:

✓ Krypton R	KRR090025	✓ Argon ²	AGN090025
✓ Krypton F	KRF090025	✓ Neon ²	NSW090025
✓ Helium	HUM090025	✓ Life F	LFE090025
✓ Xenon ² FF	XFF090025	✓ Life R	LRF090025
✓ Xenon ² H	XFF090025		
✓ Xenon ² SA	XES090025		



Back pack

Designed with special hoops to hang from the push handles or backposts of a wheelchair. Please advise users to keep their belongings secure.

Available on:

✓ Krypton R	KRR090025	✓ Life F	LFE090030
✓ Krypton F	KRF090025	✓ Life R	LRF090030
✓ Helium	HUM090030		
✓ Xenon ² FF	XFF090030		
✓ Xenon ² H	XFF090030		
✓ Xenon ² SA	XES090030		
✓ Argon ²	AGN090030		
✓ Neon ²	NSW090030		



Net back bag

Simple net bag that can be used to hold belongings at the back of the chair. Very lightweight.

Available on:

✓ Helium	HUM090020
✓ Xenon ² FF	XFF090020
✓ Xenon ² H	XFF090020
✓ Xenon ² SA	XES090020
✓ Argon ²	AGN090020
✓ Neon ²	NSW090020



ACCESSORIES

Caddy

These caddy spikes at the front of the chair can be used for carrying laptops, trays and sometimes shopping bags. Useful for users that do not have the trunk support to use a bag attached to the back of the chair.

Available on:

✓ Helium	HUM090019
✓ Xenon ² FF	XFF090019
✓ Xenon ² H	XFF090019
✓ Xenon ² SA	XES090019
✓ Argon ²	AGN090019
✓ Neon ²	NSW090019
✓ Life F	LFE090019
✓ Life R	LRF090019



Crash test option

Most chairs are fitted with brackets behind the backposts and stickers when this option is chosen. This will allow the user to be transported in the chair.

Available on:

✓ Helium	HUM090100
✓ Xenon ² FF	XFF090100
✓ Xenon ² H	XFF090100
✓ Xenon ² SA	XES090100
✓ Argon ²	AGN090100
✓ Neon ²	NSW090100
✓ Life F	LFE090100
✓ Life R	LRF090100





Mobile phone pouch

Holds the mobile at the front of the chair.

Available on:

- | | |
|-------------------------|-----------|
| ✓ Krypton R | KRR090026 |
| ✓ Krypton F | KRF090026 |
| ✓ Helium | HUM090026 |
| ✓ Xenon ² FF | XFF090026 |
| ✓ Xenon ² H | XFF090026 |
| ✓ Xenon ² SA | XES090026 |
| ✓ Argon ² | AGN090026 |
| ✓ Neon ² | NSW090026 |
| ✓ Life F | LFE090026 |
| ✓ Life R | LRF090026 |



Handbike re-enforcement tubes

For those who have opted for the handbike axle and are going to use this often we recommend the extra reinforcement tubes that attach the back of the chair to the front.

Available on:

- | | |
|-------------------------|-----------|
| ✓ Helium | HUM090014 |
| ✓ Xenon ² FF | XFF090014 |
| ✓ Xenon ² H | XFF090019 |

CASTOR FORK COLOURS

Silver Matt

Available on:

- ✓ Krypton R
- ✓ Krypton F
- ✓ Helium
- ✓ Xenon² FF
- ✓ Xenon² H
- ✓ Xenon² SA
- ✓ Neon²
- ✓ Life F
- ✓ Life R



Orange

Available on:

- ✓ Krypton R
- ✓ Krypton F
- ✓ Helium
- ✓ Xenon² FF
- ✓ Xenon² H
- ✓ Xenon² SA
- ✓ Argon²
- ✓ Life F
- ✓ Life R



Black

Available on:

- ✓ Krypton R
- ✓ Krypton F
- ✓ Helium
- ✓ Xenon² FF
- ✓ Xenon² H
- ✓ Xenon² SA
- ✓ Argon²
- ✓ Life F
- ✓ Life R



Red

Available on:

- ✓ Helium
- ✓ Xenon² FF
- ✓ Xenon² H
- ✓ Xenon² SA
- ✓ Neon²



Blue

Available on:

- ✓ Helium
- ✓ Xenon SA
- ✓ Argon²



ANODIZED WHEEL COLOURS

Blue

Available on:

- ✓ Krypton R KRR100337
- ✓ Krypton F KRF100337
- ✓ Helium HUM100337
- ✓ Xenon² FF XFF100337
- ✓ Xenon² H XFF100337
- ✓ Xenon² SA XES100337
- ✓ Neon² NSW100337



Orange

Available on:

- ✓ Krypton R KRR100336
- ✓ Krypton F KRF100336
- ✓ Helium HUM100336
- ✓ Xenon² FF XFF100336
- ✓ Xenon² H XFF100336
- ✓ Xenon² SA XES100336
- ✓ Neon² NSW100336



Red

Available on:

- ✓ Krypton R KRR100338
- ✓ Krypton F KRF100338
- ✓ Helium HUM100338
- ✓ Xenon² FF XFF100338
- ✓ Xenon² H XFF100338
- ✓ Xenon² SA XES100338
- ✓ Neon² NSW100338



ANODIZED WHEEL COLOURS

Style Package


The style package includes the option of 5 colours (Black, silver, orange, red and blue) on the castor forks, backrest bracket, axle-stem clamp and stick.

Available on:


- ✓ Krypton R
- ✓ Krypton F
- ✓ Helium
- ✓ Xenon² FF
- ✓ Xenon² H
- ✓ Xenon² SA
- ✓ Neon²



COLOUR OPTIONS



White (01)
(RAL 9010)



Signal White (32)



Sun Yellow (06)
(RAL 1018)



Tutti Fruitti
Orange (24)



Mars Red (37)



Brilliant Purple
(19)



Matt Purple (20)



Cyan Metallic (36)



Brilliant Marine
(15)



Matt Marine (16)



Brilliant Silver (02)




Matt Silver (03)



Silver (04)



Matt Anthracite
(42)



Silver Texture (26)



Matt Black (31)



Azure Blue (43)




Ruby Red (21)
(RAL 3003)




Brilliant Red (22)




Matt Red (23)



Flamingo Glitter
(41)



Pink (33)



Sky Blue (35)



Blue (14)



Kiwi (40)



Brilliant Lime
(07)




Matt Lime (08)



Graphite Metallic
(27)




Platinum Metallic
(28)



Dark Chocolate
Metallic (34)



Cappuccino (39)



Black (29)
(RAL 9005)

Note: Colours and patterns may vary based on type of chair or frame material

USEFUL RESOURCES



Eforms - <https://eforms.sunrisemedical.eu/>

- Allow you to check the latest 'no-go's'
- Show you compatible option combinations - a great training tool that develops your knowledge
- Allow you to check compatibility of options for any future customer upgrades
- Assist you in completing order forms correctly
- Generate and save customer quotations (which you can edit at a later date)
- Price your order and then enable you to save it as a PDF.

E-Sales - www.sales.sunrisemedical.co.uk

- Reserve demonstration stock online
- Check stock availability
- Buy refurbished products

Sunparts - <http://parts.sunrisemedical.eu>

- Access spare parts manuals
- Check retail prices
- Order securely online

Need a log-in? Call us on 0845 605 66 88 and we'll be happy to help!



Want to find out more about
Quickie manual chairs?

Visit
www.SunriseMedical.co.uk
where you can view the latest product
order forms, brochures and
specifications.

QUICKIE
Live without limits.



**SUNRISE
MEDICAL**



Lomax

zippie.
Explore your world!



All of our manual chairs are crash tested in accordance with ISO 7176-19.



WHITMYER

JAY
Superior Clinical Seating.



built4me

Do you have a request that's different from the options available on our regular wheelchairs?

Built-4-Me will listen to your requirements, analyse the risk factors involved and deliver a custom-made wheelchair to meet your needs.

Contact Built-4-Me:
Tel: 01384 42 15 21
Email: built4me@sunmed.co.uk



Breezy
Premium products. Quality service.

JCM
24 hour postural care



COOPERS

RGK



STERLING

SWITCH-IT
Specialist Manual & Electric Seating for People with Disabilities



Sunrise Medical Limited
Thorns Road,
Brierley Hill,
West Midlands
DY5 2LD
Tel.: +44 (0) 845 605 66 88
Fax: +44 (0) 845 605 66 89
www.SunriseMedical.co.uk

